



GUÍA DE INICIO RÁPIDO

PERSONALIZACIÓN DE LA BANDEJA DE ENTRADA DE EJEMPLO

Aprenda a:

Personaliza tu bandeja de entrada de muestras

PERSONALIZACIÓN DE LA BANDEJA DE ENTRADA DE EJEMPLO

+ Add Filter



DATE RANGE Last Year

PROGRAMA Todo

MARCA SAMPLES

OPCIONES DE IMPRESIÓN

- Select -

EXPORT XLS

INFORME(S) IMPRESO(S)



	NÚMERO LAB	GRAVEDAD	PROGRAMA	SUBPROGRAMA	CLIENTE	LUGAR DE TRABAJO	COMPONENTE	FABRICANTE DE COMPONENTE	MODELO DE COMPONENTES	SERIE DE COMPONENTES
<input type="checkbox"/>	202509222155	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	ENGINE	Heil	REL	
<input type="checkbox"/>	202509222156	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	ENGINE	Heil	REL	
<input type="checkbox"/>	202509222154	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	ENGINE	Heil	REL	
<input type="checkbox"/>	202509222153	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	ENGINE	Heil	REL	
<input type="checkbox"/>	202509222152	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	ENGINE	Heil	REL	
<input type="checkbox"/>	202509222145	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	RIGHT FRONT FINAL DRIVE	Caterpillar	730	3T30107
<input type="checkbox"/>	202509222143	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	RIGHT FRONT FINAL DRIVE	Caterpillar	730	3T30107
<input type="checkbox"/>	202509222144	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	RIGHT FRONT FINAL DRIVE	Caterpillar	730	3T30107
<input type="checkbox"/>	202509222148	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	RIGHT FRONT FINAL DRIVE	Caterpillar	730	3T30107
<input type="checkbox"/>	202509222147	B	BP CASTROL	BP CASTROL	Acme Demo	Worksite 2	RIGHT FRONT FINAL DRIVE	Caterpillar	730	3T30107

Showing 10 rows

1 of 2 << < > >>

Personaliza tu bandeja de entrada haciendo clic en el icono de engranaje y seleccionando Configuración de columna.

PERSONALIZACIÓN DE LA BANDEJA DE ENTRADA DE EJEMPLO

The screenshot displays the Castrol Labcheck interface with a 'COLUMN SETTINGS' dialog box open. The dialog has a title bar with 'COLUMN SETTINGS' and a close button. Below the title bar, there is an 'Add New Column' section with a dropdown menu currently set to '- Select -'. Underneath, the 'AVAILABLE COLUMNS' section is titled 'Drag and drop columns to reorder' and contains a list of columns, each with a drag handle icon on the left and a close 'x' icon on the right. The columns listed are: NÚMERO LAB, GRAVEDAD, PROGRAMA, SUBPROGRAMA, CLIENTE, LUGAR DE TRABAJO, COMPONENTE, FABRICANTE DE COMPONENTES, and MODELO DE COMPONENTES. A red arrow points to the 'CLIENTE' column. At the bottom of the dialog are three buttons: 'RESET', 'CANCEL', and 'SAVE'. The background shows a table with columns for 'NÚMERO LAB', 'GRAVEDAD', 'PROGRAMA', and 'SUBPROGRAMA', and a row of data with values like '202509222155', 'B', 'BP CASTROL', and 'BP CASTROL'. The top navigation bar includes 'INICIO', 'ACCIONES', 'MUESTRAS', 'ETIQUETAS', 'EQUIPO', 'ANÁLITICA', 'INFORMES', 'DISTRIBUCIÓN', 'ADMINISTRACIÓN', and 'PREFERENCIAS'. There is also a '+ Add Filter' button and a 'MARCA SAMPLES' button.

Arrastra y suelta cualquiera de las columnas para reordenar tu bandeja de entrada. Utilice la barra de desplazamiento para ver todas las columnas disponibles.

PERSONALIZACIÓN DE LA BANDEJA DE ENTRADA DE EJEMPLO

The screenshot displays the Castrol Labcheck application interface. At the top, there is a navigation menu with options: INICIO, ACCIONES, MUESTRAS, ETIQUETAS, EQUIPO, ANALÍTICA, INFORMES, DISTRIBUCIÓN, ADMINISTRACIÓN, and PREFERENCIAS. Below the menu, there is a '+ Add Filter' button and a 'DATE RANGE' dropdown set to 'Last Year'. The main data table has columns: NÚMERO LAB, GRAVEDAD, PROGRAMA, and SUBPROGRAMA. A 'COLUMN SETTINGS' dialog box is open, showing 'Add New Column' (set to '- Select -') and a list of 'AVAILABLE COLUMNS'. The available columns are: NÚMERO LAB, GRAVEDAD, PROGRAMA, SUBPROGRAMA, CLIENTE, LUGAR DE TRABAJO, COMPONENTE, FABRICANTE DE COMPONENTES, and MODELO DE COMPONENTES. A red arrow points to the 'x' icon next to 'LUGAR DE TRABAJO'. At the bottom of the dialog are 'RESET', 'CANCEL', and 'SAVE' buttons.

NÚMERO LAB	GRAVEDAD	PROGRAMA	SUBPROGRAMA
<input type="checkbox"/> 202509222155	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222156	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222154	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222153	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222152	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222145	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222143	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222144	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222148	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222147	B	BP CASTROL	BP CASTROL

Showing 10 rows

Haga clic en la x para eliminar la columna de su bandeja de entrada.

PERSONALIZACIÓN DE LA BANDEJA DE ENTRADA DE EJEMPLO

The screenshot displays the Castrol Labcheck interface. At the top, there is a navigation bar with the following menu items: INICIO, ACCIONES, MUESTRAS, ETIQUETAS, EQUIPO, ANALÍTICA, INFORMES, DISTRIBUCIÓN, ADMINISTRACIÓN, and PREFERENCIAS. Below the navigation bar, there is a '+ Add Filter' button. The main content area shows a table with columns: NÚMERO LAB, GRAVEDAD, PROGRAMA, and SUBPROGRAMA. The table contains 10 rows of sample data, all with a severity level of 'B' and a program of 'BP CASTROL'. A 'COLUMN SETTINGS' dialog box is open on the right side of the screen. The dialog has a title bar with a close button (X). Below the title bar, there is a section 'Add New Column' with a dropdown menu. The dropdown menu is open, showing a list of columns: '- Select -', '- Select -', 'Orden de trabajo' (highlighted in green), 'NÚMERO LAB', 'GRAVEDAD', 'PROGRAMA', 'SUBPROGRAMA', 'CLIENTE', 'LUGAR DE TRABAJO', 'COMPONENTE', 'FABRICANTE DE COMPONENTES', and 'MODELO DE COMPONENTES'. A red arrow points to the 'Orden de trabajo' option. At the bottom of the dialog, there are three buttons: 'RESET', 'CANCEL', and 'SAVE'. The background table shows a 'showing 10 rows' indicator at the bottom left.

NÚMERO LAB	GRAVEDAD	PROGRAMA	SUBPROGRAMA
<input type="checkbox"/> 202509222155	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222156	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222154	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222153	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222152	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222145	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222143	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222144	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222148	B	BP CASTROL	BP CASTROL
<input type="checkbox"/> 202509222147	B	BP CASTROL	BP CASTROL

Las columnas eliminadas se pueden volver a agregar haciendo clic en el menú desplegable en el campo Agregar nueva columna y seleccionando una columna.

PERSONALIZACIÓN DE LA BANDEJA DE ENTRADA DE EJEMPLO

The screenshot shows the Castrol Labcheck interface. At the top, there is a navigation bar with options: INICIO, ACCIONES, MUESTRAS, ETIQUETAS, EQUIPO, ANALÍTICA, INFORMES, DISTRIBUCIÓN, ADMINISTRACIÓN, and PREFERENCIAS. Below this is a filter section with a '+ Add Filter' button and a 'MARCA SAMPLES' button. The main area displays a table with columns: NÚMERO LAB, GRAVEDAD, PROGRAMA, and SUBPROGRAMA. The table contains 10 rows of sample data, all with a 'B' in the GRAVEDAD column. At the bottom left, it says 'Showing 10 rows'. A 'COLUMN SETTINGS' dialog box is open on the right side. It has a title bar with a close button. Below the title is an 'Add New Column' dropdown menu. Underneath is the 'AVAILABLE COLUMNS' section with the instruction 'Drag and drop columns to reorder'. The available columns are: NÚMERO LAB, GRAVEDAD, PROGRAMA, SUBPROGRAMA, CLIENTE, LUGAR DE TRABAJO, FABRICANTE DE COMPONENTES, MODELO DE COMPONENTES, and COMPONENTE. At the bottom of the dialog box are three buttons: 'RESET', 'CANCEL', and 'SAVE'. A red arrow points to the 'SAVE' button.

<input type="checkbox"/>	NÚMERO LAB	GRAVEDAD	PROGRAMA	SUBPROGRAMA
<input type="checkbox"/>	202509222155	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222156	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222154	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222153	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222152	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222145	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222143	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222144	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222148	B	BP CASTROL	BP CASTROL
<input type="checkbox"/>	202509222147	B	BP CASTROL	BP CASTROL

Para actualizar sus selecciones, seleccione el botón GUARDAR.



SERVICIO DE SOPORTE DE CASTROL LABCHECK:

Teléfono: 866-LABCHECK (522-2432)

Labchecksupport@bureauveritas.com

<https://www.labcheckresources.com/>

